

XDR Program Change Highlights

Version 3.2.14 (Jun 21st, 2018)

- Resolved issue of possible loss of exam notes when editing exams.
- Licensing service updated for better stability.

Version 3.2.13 (May 1st, 2018)

- Added ability to organize and analyze 3D scans and workup files in conjunction with Invivo and XDR's 3D Indexing service.
- Added components for launching 3D viewers.
- Added ability to import and launch 3rd party files with Windows associated applications.
- Added ability to run XDR from multiple time zones.
- Updated imaging algorithms to improve appearance of corners in radiographs.
- Added ability for XDR License Service to be accessed on computer not running the XDR database.
- Improved exam approval workflow for university installations.
- Improvements to sharpening filter to allow for future enhancements.
- Expanded range for sharpening values to allow for greater sharpening to be applied.
- Improved XDR Bridging to allow for better user feedback.
- Improved XDR Bridging to removal of bridge file after processing.
- Added ability to utilize Axium API to notify Axium of exam capture.
- Improved image display query to display images in consistent manner.
- Improved consistency of image recovery process.
- Changed temporary file location to "<user>/Appdata/roaming/XDR" unless configured differently.
- Improved XDR Licensing service to allow for FIPS Compliance.
- Updated XDR Email to support Datamotion SecureEmail API.