

# XDR Program Change Highlights

## **Version 3.2.16 (July 17<sup>th</sup>, 2020)**

- Improved image algorithms to ensure ideal diagnostic quality under less than ideal exposure
- Improved handling of Axium security group changes
- Improved handling of saving edited exams with Approver Mode activated.
- Added ability to import \*.tiff image files
- Improved notes function to ensure retention of data when copying images
- Improved image recovery to ensure images from being accidentally deleted
- Changed CAL file missing message to appear only once per exam when CAL file unavailable
- Improved aesthetic appearance of Exposure Meter
- Improved EULA message approval screen for better functionality
- Improved handling of large image files
- Updated default image treatment settings for new installations